

# Label-free Biophysical Characterization of Protein-Ion Interactions

Label-free  
MicroScale  
Thermophoresis

Thomas Schubert

2bind GmbH, Regensburg

**2bind**  
molecular interactions



# 2bind GmbH

Our mission is to give researchers, developers, and decision makers the power to accelerate and streamline key research processes. We offer tailor-made biophysics services and solutions that are perfectly adjusted to your precise needs in the key, high-value areas of Drug Discovery, Antibody Research, Aptamer Characterization, and Protein Biophysics.

A state-of-the-art portfolio of biophysical technologies enables us to provide you with a complete and comprehensive picture of your specific question: From affinity and stability to kinetics and thermodynamics. This portfolio includes MicroScale Thermophoresis (MST), nano-Differential Scanning Fluorimetry (nanoDSF), Biolayer Interferometry (BLI), and Isothermal Titration Calorimetry (ITC).

Work with the Experts in Biophysics


## Technology

MicroScale Thermophoresis (MST) is an optical fluorescent method. It records the changes in fluorescence of a target molecule as a function of temperature and the concentration of the cognate ligand molecule. Thus, in MST, the target must be fluorescent and the ligand must not be fluorescent in the same wavelength range as the target. The fluorescence changes are the result of two effects:

First, molecules exhibit directed movement along temperature gradients in solution, an effect called thermophoresis (Duhr & Braun, 2006). Such thermophoretic movement of a molecule is characterized by its size (hydrodynamic radius), its surface charge, and its hydration shell. All three parameters can be affected by binding of a ligand molecule. Thus, a concentration-dependent titration of the target molecule's ligand induces changes in the thermophoretic movement. These changes translate in quantifiable, spatial fluorescence changes, which are easily tracked optically via the target fluorescence.

Second, fluorescence is a function of the temperature. Thus, the fluorescence of the target molecule varies with the temperature (Baaske *et al.*, 2010). This temperature-dependence is additionally affected by the local molecular surroundings of the fluorophore. Thus, binding of a ligand molecule to the fluorescent target can change the chemical environment of the target fluorophore and thus the overall detected temperature-dependence of the fluorescence. This effect is known as TRIC (temperature-related intensity change; Gupta, Duhr & Baaske, 2018).

The thermophoresis and TRIC signals are additive and thus both contribute to the high sensitivity and robustness of MST measurements towards molecular binding events of all kinds. Thus, MST can be used for determining the affinity and binding strength of almost any kind of molecular interaction with very low sample consumption and very high sensitivity.



**Figure 1. Technical MST setup.** (A) MST measurements take place in small glass capillaries. Infrared and fluorescence lasers are used for generation of the MST effect and sample tracking. (B) TRIC and thermophoresis together account for a time-dependent change in fluorescence upon infrared-heating of the sample capillaries. (C) Multiple MST traces are recorded for different mixture ratios of target and ligand molecules. (D) Dose-response analysis of the MST traces allows for determination of the steady-state affinity of the target-ligand interaction.

## Features

- ✓ Wide molecule size range
- ✓ Real-time quality controls
- ✓ No immobilization required
- ✓ Very short analysis time
- ✓ Free choice of assay buffer
- ✓ Low sample consumption
- ✓ Wide affinity range
- ✓ Wide temperature range

## Advantages

- Analyze molecules from 100 Da up to 1 MDa
- Online aggregation, precipitation, and sticking controls
- Measurement is done truly in solution
- High analysis speed enables high throughput
- Biological liquids such as serum or cell lysate possible
- Minimum of only 6  $\mu$ l per sample, nM concentrations
- Affinities can be analyzed in the pM-mM range
- Analysis possible from 20°C to 45°C

## Introduction

A great number of enzymatic reactions are dependent on cofactor-binding. However, an analysis of the interaction between enzymes and their cofactors in detail is challenging. Often, it is complicated to access the basic binding parameters such as the affinity of ions using the state-of-the-art biophysical methods. Here, we demonstrate how MST can be used to test the ion-dependence of a target enzyme.

Micrococcal nuclease (MNase) is a  $\text{Ca}^{2+}$ -dependent endo-exonuclease of *Staphylococcus aureus* and a virulence factor in multiple models of infection (Kiedrowski *et al.*, 2014). Given the ion-dependence of the MNase, we developed an MST assay to test the binding of different mono- and divalent ions to MNase. As the ions do not display any fluorescence, it was possible to perform a truly label-free, in solution MST assay by using the MNase tryptophan fluorescence for MST detection ("label-free" assay).

## Results

As expected from its  $\text{Ca}^{2+}$ -dependence, MNase showed a strong preference towards divalent cations such as  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$ ,  $\text{Mn}^{2+}$ , and  $\text{Zn}^{2+}$  with  $K_D$ -values in the low millimolar concentration range ( $K_D$ :  $3.6 \pm 0.9$  mM,  $9.9 \pm 5.4$  mM,  $0.9 \pm 0.4$  mM and  $0.2 \pm 0.1$  mM). In contrast, no interaction of MNase and  $\text{Na}^+$  was detectable (Figure 2). This assay highlights not only the potential of MST to study membrane-associated proteins but also to detect binding of very small ions to large proteins; a task that is not feasible with most other biophysical methods (or requires much more effort).

Assay type	MST
Target	MNase ( <i>S. aureus</i> )
Labeling	Label-free (intrinsic protein tryptophan fluorescence).
Ligands	CaCl <sub>2</sub> ; MgCl <sub>2</sub> ; MnCl <sub>2</sub> ; ZnCl <sub>2</sub> ; NaCl
Measured K <sub>b</sub> [mM]	3.6 ± 0.9 (CaCl <sub>2</sub> ) 9.9 ± 5.4 (MgCl <sub>2</sub> ) 0.9 ± 0.4 (MnCl <sub>2</sub> ) 0.2 ± 0.1 (ZnCl <sub>2</sub> ) no binding (NaCl)




Figure 2: MST binding isotherms of MNase and various mono- and divalent ions. Data are the mean and standard deviation of two independent experiments and are fitted to a  $K_D$ -binding model under a 1:1 binding situation.

## Conclusion

MST assays are very versatile, because they can easily accommodate different buffers, detergents, as well as target and ligand types (proteins, nucleic acids, small molecules, ions). Here, we successfully established an MST assay that allows an analysis of the binding affinity of enzymes and their ionic cofactors in a fast and precise way.

Since MST features a very low sample consumption (only 6  $\mu$ l per sample is required, with concentrations of fluorescent target down to 0.5 nM), a wide molecular size range (100 Da- 1 MDa), a very short analysis time, and the strongpoint that measurements require no target immobilization, MST should be your method of choice for studying challenging targets such as membrane proteins. Moreover, MST works with virtually every available detergent, which can be important for solubilization of membrane proteins.

## Literature

Baaske *et al.* Optical thermophoresis for quantifying the buffer dependence of aptamer binding. *Angew. Chem. Int. Ed Engl.* 2010; 49(15): 2238-2241

Duhr & Braun. Why molecules move along a temperature gradient. *Proc. Natl. Acad. Sci. USA* 2006; 103(52): 19678-19682.

Gupta, Duhr & Baaske. MicroScale Thermophoresis (MST). *Encyclopedia of Biophysics*. 2018; 1-5

Kiedrowski *et al.* *Staphylococcus aureus* Nuc2 Is a Functional, Surface-Attached Extracellular Nuclease. *PLoS One*. 2014; 9(4): e95574.

Overington *et al.* How many drug targets are there? *Nat. Rev. Drug. Discov.* 2006; 5:993-996.

Rask-Andersen *et al.* Trends in the exploitation of novel drug targets. *Nat. Rev. Drug. Discov.* 2011; 10:579-590.



[2bind.com](http://2bind.com)

[info@2bind.com](mailto:info@2bind.com)

+49 941 2000 0890

Am Biopark 11  
93053 Regensburg  
Germany